



Sports Medicine

Orthopaedics & Traumatology

ALLOTRUE™ TECHNOLOGY

AlloTrue technology is a unique patented cleaning process designed to penetrate deep within donor tissue while still maintaining tissue integrity.



Technology

AlloTrue utilizes a combination of treatments of antibiotics, alcohol, peroxide (bone allografts only) and multiple water rinses. The deep penetrating cleansing, with minimal tissue reduction, results in microbial inactivation.


Process

The AlloTrue system contains motion for complete exposure of the allograft tissues to the cleansing reagents and agitation to facilitate the removal of blood and lipids. AlloTrue Cleansing Process – a unique, patented process designed to penetrate deep within donor tissue to remove blood, lipids and reduce bioburden.

PRECISION MEASURING

The length, folded diameter and single diameter are measured in every non-bone tendon. The folded diameter is measured by folding the tendon at the center over sterile umbilical tape.

The folded tendon is then passed with minimal force through sequentially smaller diameters on a sizing block until the sizing block can be slightly lifted up from the sterile field with the tendon remaining in the sizing block.



ALLOTRUE™ CLEANSING PROCESS

Created through a collaboration of AlloSource's engineering, research and development and operations departments, AlloTrue utilizes a combination of treatments of antibiotics, alcohol, peroxide (cut tissue only) and multiple water rinses. The deep cleansing, with minimal tissue manipulation, results in microbial inactivation.

- Tissue cleansing is done in a sealed, aseptic vessel utilizing an automatic process in order to maintain consistency.
- Ultrasonic cleansing reaches deep within the tissues at a tightly controlled temperature to avoid overheating tissue.
- Tissue is systematically rotated to provide maximum exposure of all tissue surfaces to cleansing reagents, avoiding the "stagnant layering" of other cleansing processes where tissue is stationary.
- Computerized controls are in place so that technicians can uphold consistent processing, rotation, temperature regulation and exposure to reagents in order to maintain tissue integrity.
- Matching the donor tissue to the reagents minimizes overexposure to cleansing agents, excessive over-processing and over-whitening.
- The microprocessor automatically determines the length of the sonication cycles and when to drain the system.
- The sonication/draining cycle is repeated several times, using new water and reagents each time.

FACT 1

AlloTrue™, a patented tissue cleansing process, removes blood and lipids and reduces bioburden from allograft tendons without exposing the tissue to harsh chemicals.

FACT 2

Controlled low dose (10-15 KGY), low temperature irradiation does not impact biomechanical properties of allograft tendon.

Tissues maintain structural and biomechanical properties similar to those of non-irradiated tendon allograft.

Allosource uses a validated low irradiation dose (< 10-15) kilogray and low temperature which preserves biomechanical properties of tissue.

FACT 3

Controlled low dose (10-15 KGY), low temperature irradiation does not impact clinical results of allograft tendons.

Numerous clinical studies support the use of controlled, low dose low temperature irradiation of allograft tendons, and provide optimal clinical results of sterile allografts.



Semitendon

AS44317003	Semitendon (Single Strand)	21–23 cm
1824-14	Semitendon (Single Strand)	21–23 cm
1824-14XL*	Semitendon (Single Strand)	>24 cm
1014-14	Semitendon	21–23 cm
1014-14XL*	Semitendon	>24 cm

Peroneus Longus

AS43917000	Peroneus Longus	26–30 cm
1170-14	Peroneus Longus	26–30 cm
AS43917000XL*	Peroneus Longus	>30 cm
1170-14XL*	Peroneus Longus	>30 cm
1170-14S*	Peroneus Longus	<23 cm

Gracilis

AS44317002	Gracilis (Single Strand)	18–23 cm
3167-14	Gracilis (Single Strand)	18–23 cm
1009-14XL*	Gracilis	>25 cm



Tibialis Anterior

AS44317000	Tibialis Anterior (Single Strand)	21–23 cm
AS41517000	Tibialis Anterior (Single Strand)	>23 cm

Tibialis Posterior

AS44317001	Tibialis Posterior (Single Strand)	21–23 cm
AS41617000	Tibialis Posterior (Single Strand)	>23 cm




Achilles WITH bone block

AS10017000	Achilles
1025-14	Achilles
1025-14XL*	Achilles XL

>24 cm

Achilles WITHOUT bone block

1085-14	Achilles
---------	----------



Quadriceps

1903-14	Quadriceps with Bone
---------	----------------------



Femoral Head frozen

AS14817048	Femoral Head 1/1
1005-14	Femoral Head 1/1
NT31110103	Femoral Head 1/2



Patella Ligament with Extensor

1162-14R

Patella Lgmt. with Extensor Right

1162-14L

Patella Lgmt. with Extensor Left



Patella Ligament Whole

AS17917000

Patella Lgmt. Whole



Patella Ligament Hemi

AS17817000

Patella Lgmt. Hemi

>4,2 cm



Patella Ligament Pre-shaped

1964-14

Patella Lgmt. Pre-shaped 1,0 cm

3398-14

Patella Lgmt. Pre-shaped 1,2 cm

PROCEDURE FOR ORDERING A MENISCUS ALLOGRAFT


As a first step we ask you to send us a CD (in duplicate) with the patient's MRI images.

These data will be forwarded to our supplier / bone bank. Here the images will be measured in order to find a suitable meniscus.

There is no definable delivery time. It depends on the availability of the measured article.

If a suitable meniscus allograft is found, the treating physician receives a document with the details of the graft. This document should be reviewed and signed to serve as an order confirmation. It is important that a cost approval from the insurance company is already available at this time. The meniscus cannot be returned after signature.

Now the order is placed with our supplier and the graft is reserved exclusively. The delivery to the hospital is made on the surgery date by Innight in a box with dry ice.



Meniscus

AS28325001 Lateral Right
90B011-X Lateral Right

AS28325002 Lateral Left
90B008-X Lateral Left

AS28225001 Medial Right
90B010-X Medial Right

AS28225002 Medial Left
90B007-X Medial Left



ORIGIN OF THE TRANSPLANTS

USA / EU



STORAGE

Musculoskeletal tissue may be stored at or below -40°C.



SHIPPING

The allografts are delivered in a carton filled with dry ice by Innight or Swiss Med Services.



APPLICATION

To thaw, place the innermost pouch in a sterile basin on a sterile field. Immerse inner pouch completely in warm, sterile isotonic solution. The recommended thawing time is 30–60 minutes, depending upon the size of the graft. If removal of blood and marrow elements is desired, rinse graft completely using high-pressure lavage with an isotonic solution.

SWISS MEDTECH

Member

Neutromedics AG is a member of Swiss Medtech.



Neutromedics AG is ISO 13485 certified.



Schweizerische Eidgenossenschaft
Confédération suisse
Confederazione Svizzera
Confederaziun svizra

BAG permit: Import & Export of Human Tissues

NEUTROMEDICS

NEUTROMEDICS AG
Alte Steinhauserstrasse 19
CH-6330 Cham/Switzerland

Tel. +41 (0)41 748 04 40
Fax +41 (0)41 748 04 41
e-mail: info@neutromedics.ch
www.neutromedics.ch



Klimaneutral
Druckprodukt
ClimatePartner.com/15421-2011-1001

