

BATH COMMODE CHAIR

moem

*Experts in well being
of people*




www.tecnimoem.com


tecnimoem
Research+Development+Innovation



MOEM bath commode chair



With MOEM bath commode chair will reduce people's movements facilitating the cleanliness



WEIGHT
15 kg



MAXIMUM
USER WEIGHT
150 kg

Materials


- » Steel frame painted with epoxy polyester
- » Plastic parts with polypropylene and ABS
- » Stainless steel screws
- » Poliamide wheels

FEATURES


- » Removable armrest. Blockable.
- » Light and easy to manoeuvre.
- » Removable cover.
- » Opening at the back to facilitate access and clean.
- » Footrest, five position.
- » Ergonomic handle.
- » Ergonomic back, possibility to add seat belt.
- » Epoxy polyester paint.

- » Prevent liquid penetration.
- » Individual brake.
- » 360° rotation.
- » Easy assembly, without tools.
- » Adaptable in any toilet.
- » Removable chamber pot.
- » Make easy the cleanliness of person with disabilities.


TECHNICAL DETAILS



*Ergonomic
back*



*Ergonomic handle
Easy to manouevre*



*Removable cover
More comfort for
user.
Facilitates access
and clean.*



*Removable armrest
Blockable*



*Individual brake
On the four wheels*



*Footrest
Five position.
Fully collapsible*





www.tecnimoem.com

Tecnimoem 97, S.L. Ctra. Logroño-Mendavia,
Km. 5 - 31230 Viana (Navarra)
Tel. +34 948 646 213 - Fax. +34 948 646 380
tecnimoem@tecnimoem.com



Management
System
ISO 9001:2008
www.tuv.com
ID 910506418

Licencia de fabricación
de Producto Sanitario
Nº 5867-PS

